




96-well microplates BRANDplates® with transparent bottom

Polystyrene, untreated or optimized for immunological and cell culture applications.

Cat.No.	Bottom	Surface	Color	Type	Sterile	max. Vol.	Lid
7816 10	F	pureGrade™	white	transp. bottom	-	330 µl	-
7816 11	F	pureGrade™	black	transp. bottom	-	330 µl	-
7816 70	F	pureGrade™ S	white	transp. bottom	+	330 µl	+
7816 71	F	pureGrade™ S	black	transp. bottom	+	330 µl	+
7817 31	F	immunoGrade™	black	transp. bottom	-	330 µl	-
7818 50	F	lipoGrade™	white	transp. bottom	-	330 µl	-
7818 51	F	lipoGrade™	black	transp. bottom	-	330 µl	-
7819 70	F	cellGrade™	white	transp. bottom	+	330 µl	+
7819 71	F	cellGrade™	black	transp. bottom	+	330 µl	+
7820 30	F	cellGrade™ plus	white	transp. bottom	+	330 µl	+
7820 31	F	cellGrade™ plus	black	transp. bottom	+	330 µl	+
7820 90	F	cellGrade™ premium	white	transp. bottom	+	330 µl	+
7820 91	F	cellGrade™ premium	black	transp. bottom	+	330 µl	+
7819 10	F	inertGrade™	white	transp. bottom	+	330 µl	+
7819 11	F	inertGrade™	black	transp. bottom	+	330 µl	+

untreated
immunological
cell culture

All BRANDplates® microplates are manufactured according to the requirements of the ANSI-/SLAS-standards 1 to 4 - 2004. This ensures the compatibility of the BRANDplates® when performing automated processes. Subject to technical modification without notice. Errors excepted.